EASIROC User's Guide

December 10, 2019

1 About EASIROC

We can construct a tightly integrated detector with monolithic MPPC array. However, to cover the large area of the detector, we need to handle large amount of MPPCs. MPPCs on the single monolithic MPPS array have different break down voltage and gain, thus it is necessary to adjust the operation voltage for each MPPC.

To much these needs, KEK and openIt developed EASIROC board (GN-1101-1,GN-1101-2R) in 2013. This chip has some characteristics listed below:

- Readout 32 channels simultaneously
- Take positive voltage as input, amplifying, shaping output positive voltage
- Adjust bias voltage by 8 bit (0-4.5 V)
- Have discriminator in it
- Set parameters through slow control
- Cover the dynamics of 160 fC - 320 pC as input charge (this corresponds to 1-2000 p.e. with the assumption that MPPC's gain is 10^6)

University of Tokyo upgraded the EASIROC for the use of MPPC mass test for WAGASCI experiment and attached the function for scaler and TDC information.

2 Internal Operation of EASIROC

The board has 2 EASIROC chips and can read out 64 channels at the same time at maximum. The analog signals from EASIROC chips are converted to digital signal with 4 ADC. Discriminator output is sent to MHTDC and Scaler, then transported to data taking PC (DAC) through Gatherer, Sender, SiTCP.

DAC and EASIROC board are connected with Ethernet cable. The board's operation voltage is +6 V and it can be supplied from NIM crate.


Figure 1: Side-view of NIM-EASIROC board


Figure 2: Circuit diagram of NIM-EASIROC board


Figure 3: Diagram of EASIROC chip

The way of converting the analog signal to the digital one is like this:

1. Integrate and amplify the charge by pre-amplifier which has two types of gain (low and high).
2. Send the pre-amplifier output to fast-shaper and slow-shaper.
3. Signal sent to the fast-shaper is shaped with the time constant 15 ns and converted to the digital signal if it is over the threshold of discriminator.
4. Signal sent to the slow-shaper is shaped with the time constant 25-175 ns and recorded its voltage value when it get HOLD signal.

3 How to Use EASIROC

3.1 Trigger signals

Basically, only external trigger can be used. When you want to use self trigger, use output signals with appropriate delay. In that case, you should enter the signal with the order; "HOLD", "T STOP", "ACCEPT". You should insert more than 2 us between HOLD and ACCEPT. It depends on ADC converting rate but 2.5 us is maybe enough.


Figure 4: An example of self trigger circuit[1]

Figure 5: The front panel of EASIROC board

3.2 Interactive mode

EASIROC has an inner read line so you can handle EASIROC module interactively.

3.2.1 Handle interactive mode

Start interactive mode.

```
$ ./Controller.rb [IP Address (default: 192.168.10.16)]
```

Stop interactive mode.

```
> exit [quit]
```

3.2.2 Set bias voltage

Print the voltage and current value of HV.

```
> statusHV
```

Set the HV at the value of [bias voltage].

```
> setHV [bias voltage]
```

Increase HV with some steps. In each step, check the current and stop if it reaches to the limit.

```
> increaseHV [bias voltage]
```

3.2.3 Data taking

Reset the slow control values.

```
> slowcontrol
```

Start data acquisition.

```
> read [Event #] [Filename]
```

ON/OFF the ADC [TDC/Scaler].

```
> adc [on/off]
> tdc [on/off]
> scaler [on/off]
```

It makes data taking faster if you off the ADC [TDC/Scaler] when you don't need them.

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3.3.1 RegisterValue.yaml

```
EASIROC1: # Set for each chip
  Capacitor HG PA Fdbck: 100fF # Define amplifying level. Capacitor values
  Capacitor LG PA Fdbck: 100fF # are contained in RegisterValueAlias.yaml
  Time Constant HG Shaper: 100ns # Time constant of slow-shaper
  Time Constant LG Shaper: 50ns
  DAC code: 600 # Threshold of the discriminator after fast-shaper

EASIROC2:
  Capacitor HG PA Fdbck: same
  Capacitor LG PA Fdbck: same
  Time Constant HG Shaper: same
  Time Constant LG Shaper: same
  DAC code: same

High Gain Channel 1: 0 # Read out channel from HG1/HG2
High Gain Channel 2: -1 # 0 : read out || 1 : not read out
Probe Channel 1: -1 # Output from the probe of front panel
Probe Channel 2: -1
Probe 1: Out_fs
Probe 2: Out_fs # Out_PA_HG, Out_PA_LG, Out_ssh_HG, Out_ssh_LG, Out_fs
SelectableLogic:
  Pattern: Or64 # OneCh_#, Or32u, Or32d, Or64, Or32And,...
  HitNum Threshold: 4 # Threshold for each OR logic. 0~64. Default: 0
  And Channels: -1 # Cannels used in And Logic. 0~63. Default: -1
TimeWindow: 4095ns
UsrClkOut: "OFF" # Periodic signal from syn out of front panel # "ON", 1Hz,...
Trigger: ## This "Trigger" values are not used for this version.
  Mode: 0 #0-7
  DelayTrigger: -1 #500MHz #default:-1, 0-253 #trig -> hold -> l1 -> l2
  DelayHold: -1 #25MHz
  DelayL1Trig: -1 #6MHz
  Width: raw
```

When increase the DAC of discriminator, threshold get lower.

3.3.2 InputDAC.yaml

There written the 8-bit Input DAC values for 32 channels by chip. They can be changed between 256-511 (The top bit is always set to 1 (=enable)).

When increase the DAC value, bias voltage get lower.

```
---
EASIROC1:
  Input 8-bit DAC:
    - 350
    - 350
    - 350
    - 350
    - 350
    - 350
    - 350
    - 350
    - 350
```

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